

**No. 6400-00 HANDY TYPE DIGITAL PH METER MODEL SK-620PH****SPECIFICATIONS**

Measuring range:	pH: 0.00 to 14.00pH Temp.: 0.0 to 50.0°C
Accuracy:	pH: $\pm(0.3+1 \text{ digit})\text{pH}$ at 0 to 2pH $\pm(0.1+1 \text{ digit})\text{pH}$ at 2 to 12pH $\pm(0.4+1 \text{ digit})\text{pH}$ at 12 to 14pH Temp.: $\pm 0.5^\circ\text{C}$ at 20 to 30°C, $\pm 0.7^\circ\text{C}$ at other
Resolution:	0.01pH / 0.1°C
Elements:	pH: Glass electrode Temp.: Thermistor
Power requirements:	9 VDC (6F22 (006P) x 1 pc.) or AC adapter (optional)
Materials:	Case: ABS resin Probe: Glass
Dimensions:	72(W) x 182(H) x 31.5(D) mm
Weight:	about 220 g (including battery)
Standard accessories:	PHP-31 probe, bottle for glass electrode, vinyl cover, carrying case, 006P battery x 1, Standard liquid set (Phthalate solution, Neutral phosphate solution, Carbonate solution)

**No. 6410-00 PENTYPE DIGITAL PH METER MODEL SK-610PHII****SPECIFICATIONS**

Measuring range:	pH: 2.0 to 12.0pH Temp.: 0.0 to 50.0°C
Accuracy:	pH: $\pm 2\text{pH}$ Temp.: $\pm 0.6^\circ\text{C}$ at 20 to 30°C, $\pm 0.8^\circ\text{C}$ at other
Resolution:	0.01pH / 0.1°C
Elements:	pH: Glass electrode Temp.: Thermistor
Power requirements:	6VDC (3V Lithium Battery CR2032 x 2)
Materials:	Case: ABS resin Probe: Glass
Dimensions:	37(W) x 150(H) x 42(D) mm
Weight:	about 77 g (including batteries)
Standard accessories:	Vinyl pouch, lithium battery x 2

**POCKET TYPE DIGITAL PH METERS****No. 6416-00 MODEL SK-630PH****No. 6418-00 MODEL SK-631PH****No. 6420-00 MODEL SK-632PH****FEATURES**

- Calibration function
- HOLD function
- Auto power-off function
- Auto temperature compensation  
(Not applicable for model SK-630PH)



SK-630PH



SK-631PH



SK-632PH

**SPECIFICATIONS**

Model		SK-630PH	SK-631PH	SK-632PH
Measuring range	pH	2.0 to 12.0pH	2.0 to 12.0pH	2.00 to 12.00pH
	Temp.	—	0.0 to 50.0°C	0.0 to 50.0°C
Accuracy	pH	±0.6pH	±0.5pH	±0.4pH
	Temp.	—	±1.0°C	±0.5°C
Resolution	pH	0.1pH	0.1pH	0.01pH
	Temp.	—	0.5°C	0.1°C
Elements	pH	Glass electrode		
	Temp.	—	Thermistor	
Measuring ambient		0 to 50°C		
Power requirement		6VDC (Alkaline Battery LR44 1.5VDC x 4)		
Battery life		about 20 hours in continuous use		
Materials		Case: ABS resin    Sensor: Glass		
Dimensions		36 (W) x 150 (H) x 25 (D) mm		
Weight		about 65 g (including batteries)		
Standard accessories		Alkaline Battery LR44 (GPS76E) x 4		

